

Human ATXN10/Ataxin-10 ELISA Kit

Catalog No: #EK19716

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ATXN10/Ataxin-10 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Ataxin-10
Target Name	Ataxin-10
Other Names	Brain protein E46 homolog, Spinocerebellar ataxia type 10 protein, SCA10
Accession No.	Q9UBB4
Uniprot	Q9UBB4
GeneID	25814;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.8pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ATXN10

Note: This product is for in vitro research use only