

Rat Amhr2/Anti-Muellerian hormone type-2 receptor ELISA Kit



Catalog No: #EK19717

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Amhr2/Anti-Muellerian hormone type-2 receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Anti-Muellerian hormone type-2 receptor
Target Name	Anti-Muellerian hormone type-2 receptor
Other Names	Anti-Muellerian hormone type II receptor,AMH type II receptor,C14,MIS type II receptor,MISR11
Accession No.	Q62893
Uniprot	Q62893
GeneID	29530;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:30pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Amhr2

Note: This product is for in vitro research use only