

# Mouse Crem/cAMP-responsive element modulator ELISA Kit

Catalog No: #EK19718

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Crem/cAMP-responsive element modulator ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse cAMP-responsive element modulator
Target Name	cAMP-responsive element modulator
Other Names	Inducible cAMP early repressor, ICER
Accession No.	P27699
Uniprot	P27699
GeneID	12916;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Crem

Note: This product is for in vitro research use only