

## Human CTRB1/Chymotrypsinogen B ELISA Kit

Catalog No: #EK19719

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CTRB1/Chymotrypsinogen B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Chymotrypsinogen B
Target Name	Chymotrypsinogen B
Other Names	CTRB
Accession No.	P17538
Uniprot	P17538
GeneID	1504;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.4pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CTRB1

Note: This product is for in vitro research use only