

Rat Ctrb1/Chymotrypsinogen B ELISA Kit

Catalog No: #EK19720

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Ctrb1/Chymotrypsinogen B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Chymotrypsinogen B
Target Name	Chymotrypsinogen B
Other Names	Ctrb
Accession No.	P07338
Uniprot	P07338
GeneID	24291;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.055ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ctrb1

Note: This product is for in vitro research use only