

## Human BTB/Biotinidase ELISA Kit

Catalog No: #EK19721

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human BTB/Biotinidase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Biotinidase
Target Name	Biotinidase
Other Names	Biotinase
Accession No.	P43251
Uniprot	P43251
GeneID	686;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 18pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: BTB

Note: This product is for in vitro research use only