

Human CD8B/T-cell surface glycoprotein CD8 beta chain ELISA Kit

Catalog No: #EK19722

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CD8B/T-cell surface glycoprotein CD8 beta chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human T-cell surface glycoprotein CD8 beta chain
Target Name	T-cell surface glycoprotein CD8 beta chain
Other Names	CD8b,CD8B1
Accession No.	P10966
Uniprot	P10966
GeneID	926;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 24pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CD8B

Note: This product is for in vitro research use only