

## Human ADK/Adenosine kinase ELISA Kit

Catalog No: #EK19723

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ADK/Adenosine kinase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Adenosine kinase
Target Name	Adenosine kinase
Other Names	AK, Adenosine 5'-phosphotransferase
Accession No.	P55263
Uniprot	P55263
GeneID	132;

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.3ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: ADK

Note: This product is for in vitro research use only