

Human CNR1/Cannabinoid receptor 1 ELISA Kit

Catalog No: #EK19724

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CNR1/Cannabinoid receptor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cannabinoid receptor 1
Target Name	Cannabinoid receptor 1
Other Names	CB-R,CANN6,CNR
Accession No.	P21554
Uniprot	P21554
GeneID	1268;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.092ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CNR1

Note: This product is for in vitro research use only