

Rat Cnr1/Cannabinoid receptor 1 ELISA Kit

Catalog No: #EK19725

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cnr1/Cannabinoid receptor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cannabinoid receptor 1
Target Name	Cannabinoid receptor 1
Other Names	CB-R,Brain-type cannabinoid receptor,Skr6
Accession No.	P20272
Uniprot	P20272
GeneID	25248;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.089ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cnr1

Note: This product is for in vitro research use only