

# Human BTLA/B- and T-lymphocyte attenuator ELISA Kit



Catalog No: #EK19726

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human BTLA/B- and T-lymphocyte attenuator ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human B- and T-lymphocyte attenuator
Target Name	B- and T-lymphocyte attenuator
Other Names	B- and T-lymphocyte-associated protein,CD272
Accession No.	Q7Z6A9
Uniprot	Q7Z6A9
GeneID	151888;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 19pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: BTLA

Note: This product is for in vitro research use only