

Mouse Adnp/Activity-dependent neuroprotector homeobox protein ELISA Kit

Catalog No: #EK19727

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Adnp/Activity-dependent neuroprotector homeobox protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Activity-dependent neuroprotector homeobox protein
Target Name	Activity-dependent neuroprotector homeobox protein
Other Names	Activity-dependent neuroprotective protein
Accession No.	Q9Z103
Uniprot	Q9Z103
GeneID	11538;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.053ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adnp

Note: This product is for in vitro research use only