

# Human BNIP3/BCL2/adenovirus E1B 19 kDa protein-interacting protein 3 ELISA Kit

Catalog No: #EK19728

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human BNIP3/BCL2/adenovirus E1B 19 kDa protein-interacting protein 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human BCL2/adenovirus E1B 19 kDa protein-interacting protein 3
Target Name	BCL2/adenovirus E1B 19 kDa protein-interacting protein 3
Other Names	NIP3
Accession No.	Q12983
Uniprot	Q12983
GeneID	664;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 18pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: BNIP3

Note: This product is for in vitro research use only