

# Mouse Bnip3/BCL2/adenovirus E1B 19 kDa protein-interacting protein 3 ELISA Kit

Catalog No: #EK19729

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Bnip3/BCL2/adenovirus E1B 19 kDa protein-interacting protein 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse BCL2/adenovirus E1B 19 kDa protein-interacting protein 3
Target Name	BCL2/adenovirus E1B 19 kDa protein-interacting protein 3
Other Names	Nip3
Accession No.	O55003
Uniprot	O55003
GeneID	12176;

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.055ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Bnip3

Note: This product is for in vitro research use only