

Human AKR7A2/Aflatoxin B1 aldehyde reductase member 2 ELISA Kit



Catalog No: #EK19730

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human AKR7A2/Aflatoxin B1 aldehyde reductase member 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Aflatoxin B1 aldehyde reductase member 2
Target Name	Aflatoxin B1 aldehyde reductase member 2
Other Names	AFB1 aldehyde reductase 1, AFB1-AR 1, Aldoketoreductase 7, Succinic semialdehyde reductase, SSA reductase, AFAR, AFAR1, AKR7
Accession No.	O43488
Uniprot	O43488
GeneID	8574;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.069ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: AKR7A2

Note: This product is for in vitro research use only