

Mouse Akr7a2/Aflatoxin B1 aldehyde reductase member 2 ELISA Kit

Catalog No: #EK19731

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Akr7a2/Aflatoxin B1 aldehyde reductase member 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Aflatoxin B1 aldehyde reductase member 2
Target Name	Aflatoxin B1 aldehyde reductase member 2
Other Names	Succinic semialdehyde reductase,SSA reductase,Afar,Akr7a5
Accession No.	Q8CG76
Uniprot	Q8CG76
GeneID	110198;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Akr7a2

Note: This product is for in vitro research use only