

Rat Akr7a2/Aflatoxin B1 aldehyde reductase member 2 ELISA Kit



Catalog No: #EK19732

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Akr7a2/Aflatoxin B1 aldehyde reductase member 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Aflatoxin B1 aldehyde reductase member 2
Target Name	Aflatoxin B1 aldehyde reductase member 2
Other Names	rAFAR2,Succinic semialdehyde reductase,SSA reductase,Afar2,Aiar
Accession No.	Q8CG45
Uniprot	Q8CG45
GeneID	171445;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.17ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Akr7a2

Note: This product is for in vitro research use only