

RALGDS Antibody

Catalog No: #44681

Package Size: #44681-1 50ul #44681-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RALGDS Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	RALGDS antibody detects endogenous levels of total RALGDS
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human RALGDS
Target Name	RALGDS
Other Names	ral guanine nucleotide dissociation stimulator antibody RalGEF antibody RGF antibody
Accession No.	Swiss-Prot#: Q12967 NCBI Gene ID:5900
Uniprot	Q12967
GeneID	5900;
Calculated MW	100 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Images



Western Blot analysis of Hela lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only