

Human BHMT2/S-methylmethionine--homocysteine S-methyltransferase BHMT2 ELISA Kit



Catalog No: #EK19733

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human BHMT2/S-methylmethionine--homocysteine S-methyltransferase BHMT2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human S-methylmethionine--homocysteine S-methyltransferase BHMT2
Target Name	S-methylmethionine--homocysteine S-methyltransferase BHMT2
Other Names	SMM-hcy methyltransferase, Betaine--homocysteine S-methyltransferase 2
Accession No.	Q9H2M3
Uniprot	Q9H2M3
GeneID	23743;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.081ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: BHMT2

Note: This product is for in vitro research use only