

Mouse Cfl1/Cofilin-1 ELISA Kit

Catalog No: #EK19734

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cfl1/Cofilin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cofilin-1
Target Name	Cofilin-1
Other Names	Cofilin, non-muscle isoform
Accession No.	P18760
Uniprot	P18760
GeneID	12631;

Application Details

Assay Type: Sandwich Detect Range: 6.25-400pg/mL Sensitivity: 3pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cfl1

Note: This product is for in vitro research use only