

## Human CREBRF/CREB3 regulatory factor ELISA Kit

Catalog No: #EK19735

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CREBRF/CREB3 regulatory factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human CREB3 regulatory factor
Target Name	CREB3 regulatory factor
Other Names	Luman recruitment factor,LRF,C5orf41
Accession No.	Q8IUR6
Uniprot	Q8IUR6
GeneID	153222;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.164ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CREBRF

Note: This product is for in vitro research use only