

Mouse Adamts20/A disintegrin and metalloproteinase with thrombospondin motifs 20 ELISA Kit



Catalog No: #EK19736

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Adamts20/A disintegrin and metalloproteinase with thrombospondin motifs 20 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse A disintegrin and metalloproteinase with thrombospondin motifs 20
Target Name	A disintegrin and metalloproteinase with thrombospondin motifs 20
Other Names	ADAM-TS 20
Accession No.	P59511
Uniprot	P59511
GeneID	223838;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adamts20

Note: This product is for in vitro research use only