

Human CCR1/C-C chemokine receptor type 1 ELISA Kit



Catalog No: #EK19737

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CCR1/C-C chemokine receptor type 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human C-C chemokine receptor type 1
Target Name	C-C chemokine receptor type 1
Other Names	C-C CKR-1, HM145, LD78 receptor, Macrophage inflammatory protein 1-alpha receptor, MIP-1alpha-R, RANTES-R, CD191, CMKBR1, CMKR1, SCYAR1
Accession No.	P32246
Uniprot	P32246
GeneID	1230;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.398ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCR1

Note: This product is for in vitro research use only