

Mouse Ccr1/C-C chemokine receptor type 1 ELISA Kit

Catalog No: #EK19738

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ccr1/C-C chemokine receptor type 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse C-C chemokine receptor type 1
Target Name	C-C chemokine receptor type 1
Other Names	C-C CKR-1, Macrophage inflammatory protein 1-alpha receptor, MIP-1alpha-R, RANTES-R, CD191, Cmkbr1
Accession No.	P51675
Uniprot	P51675
GeneID	12768;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/ml Sensitivity: 0.399ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ccr1

Note: This product is for in vitro research use only