

Human AQP12A/Aquaporin-12A ELISA Kit

Catalog No: #EK19739

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human AQP12A/Aquaporin-12A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Aquaporin-12A
Target Name	Aquaporin-12A
Other Names	AQP-12,AQP12,AQPX2
Accession No.	Q8IXF9
Uniprot	Q8IXF9
GeneID	375318;

Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/ml Sensitivity:0.36ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AQP12A

Note: This product is for in vitro research use only