

## Mouse Cth/Cystathionine gamma-lyase ELISA Kit

Catalog No: #EK19740

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cth/Cystathionine gamma-lyase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cystathionine gamma-lyase
Target Name	Cystathionine gamma-lyase
Other Names	Cysteine-protein sulfhydrase,Gamma-cystathionase
Accession No.	Q8VCN5
Uniprot	Q8VCN5
GeneID	107869;

## Application Details

Assay Type: Sandwich Detect Range:6.25-400ng/mL Sensitivity:3.0ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Cth

Note: This product is for in vitro research use only