

Bovine APCS/Serum amyloid P-component ELISA Kit

Catalog No: #EK19742

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine APCS/Serum amyloid P-component ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Serum amyloid P-component
Target Name	Serum amyloid P-component
Other Names	SAP
Accession No.	Q3T004
Uniprot	Q3T004
GeneID	514488;

Application Details

Assay Type: Sandwich Detect Range: 6.25-400ng/mL Sensitivity: 3.175ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: APCS

Note: This product is for in vitro research use only