

## Human APCS/Serum amyloid P-component ELISA Kit

Catalog No: #EK19743

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human APCS/Serum amyloid P-component ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Serum amyloid P-component
Target Name	Serum amyloid P-component
Other Names	SAP,9.5S alpha-1-glycoprotein,PTX2
Accession No.	P02743
Uniprot	P02743
GeneID	325;

## Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:8.9pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:APCS

Note: This product is for in vitro research use only