

Mouse Apcs/Serum amyloid P-component ELISA Kit

Catalog No: #EK19744

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Apcs/Serum amyloid P-component ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Serum amyloid P-component
Target Name	Serum amyloid P-component
Other Names	SAP,Ptx2,Sap
Accession No.	P12246
Uniprot	P12246
GeneID	20219;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:40.6pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Apcs

Note: This product is for in vitro research use only