

Rat Apcs/Serum amyloid P-component ELISA Kit

Catalog No: #EK19745

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Apcs/Serum amyloid P-component ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Serum amyloid P-component
Target Name	Serum amyloid P-component
Other Names	SAP,Ptx2,Sap
Accession No.	P23680
Uniprot	P23680
GeneID	29339;

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 3.7pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Apcs

Note: This product is for in vitro research use only