

# Human ARHGAP21/Rho GTPase-activating protein 21 ELISA Kit



Catalog No: #EK19746

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ARHGAP21/Rho GTPase-activating protein 21 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Rho GTPase-activating protein 21
Target Name	Rho GTPase-activating protein 21
Other Names	Rho GTPase-activating protein 10,Rho-type GTPase-activating protein 21,ARHGAP10,KIAA1424
Accession No.	Q5T5U3
Uniprot	Q5T5U3
GeneID	57584;

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:34pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ARHGAP21

Note: This product is for in vitro research use only