

Mouse Arhgap21/Rho GTPase-activating protein 21 ELISA Kit



Catalog No: #EK19747

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Arhgap21/Rho GTPase-activating protein 21 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Rho GTPase-activating protein 21
Target Name	Rho GTPase-activating protein 21
Other Names	Rho GTPase-activating protein 10,Rho-type GTPase-activating protein 21,Arhgap10,Kiaa1424
Accession No.	Q6DFV3
Uniprot	Q6DFV3
GeneID	71435;

Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:10pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Arhgap21

Note: This product is for in vitro research use only