

## Dog CCL17/C-C motif chemokine 17 ELISA Kit

Catalog No: #EK19748

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dog CCL17/C-C motif chemokine 17 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog C-C motif chemokine 17
Target Name	C-C motif chemokine 17
Other Names	CC chemokine TARC, Small-inducible cytokine A17, Thymus and activation-regulated chemokine, TARC
Accession No.	Q95N01
Uniprot	Q95N01
GeneID	403586;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 40.7pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CCL17

Note: This product is for in vitro research use only