

Human CCL17/C-C motif chemokine 17 ELISA Kit

Catalog No: #EK19749

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CCL17/C-C motif chemokine 17 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human C-C motif chemokine 17
Target Name	C-C motif chemokine 17
Other Names	CC chemokine TARC, Small-inducible cytokine A17, Thymus and activation-regulated chemokine, SCYA17, TARC
Accession No.	Q92583
Uniprot	Q92583
GeneID	6361;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.093ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCL17

Note: This product is for in vitro research use only