

Human CCL23/C-C motif chemokine 23 ELISA Kit

Catalog No: #EK19750

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CCL23/C-C motif chemokine 23 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human C-C motif chemokine 23
Target Name	C-C motif chemokine 23
Other Names	CK-beta-8,CKB-8,Macrophage inflammatory protein 3,MIP-3,Myeloid progenitor inhibitory factor 1,MPIF-1,Small-inducible cytokine A23,MIP3,MPIF1,SCYA23
Accession No.	P55773
Uniprot	P55773
GeneID	6368;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/ML Sensitivity:15pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CCL23

Note: This product is for in vitro research use only