

## Human CCL21/C-C motif chemokine 21 ELISA Kit

Catalog No: #EK19751

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CCL21/C-C motif chemokine 21 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human C-C motif chemokine 21
Target Name	C-C motif chemokine 21
Other Names	6Ckine,Beta-chemokine exodus-2,Secondary lymphoid-tissue chemokine,SLC,Small-inducible cytokine A21,SCYA21,UNQ784/PRO1600
Accession No.	O00585
Uniprot	O00585
GeneID	6366;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.182ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CCL21

Note: This product is for in vitro research use only