

Mouse Ccl21a/C-C motif chemokine 21a ELISA Kit

Catalog No: #EK19752

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ccl21a/C-C motif chemokine 21a ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse C-C motif chemokine 21a
Target Name	C-C motif chemokine 21a
Other Names	6Ckine,Beta-chemokine exodus-2,Small-inducible cytokine A21a,Thymus-derived chemotactic agent 4,TCA4,Scya21,Scya21a
Accession No.	P84444
Uniprot	P84444
GeneID	100504362;18829;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.184ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ccl21a

Note: This product is for in vitro research use only