

Rat Cd2/T-cell surface antigen CD2 ELISA Kit

Catalog No: #EK19753

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cd2/T-cell surface antigen CD2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat T-cell surface antigen CD2
Target Name	T-cell surface antigen CD2
Other Names	LFA-2,LFA-3 receptor,OX-34 antigen,T-cell surface antigen T11/Leu-5,CD2
Accession No.	P08921
Uniprot	P08921
GeneID	497761;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.186ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name: Cd2

Note: This product is for in vitro research use only