

Human ARHGEF1/Rho guanine nucleotide exchange factor 1 ELISA Kit



Catalog No: #EK19754

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Human ARHGEF1/Rho guanine nucleotide exchange factor 1 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Hu |
| Specificity | Natural and recombinant Human Rho guanine nucleotide exchange factor 1 |
| Target Name | Rho guanine nucleotide exchange factor 1 |
| Other Names | 115 kDa guanine nucleotide exchange factor,p115-RhoGEF,Sub1.5 |
| Accession No. | Q92888 |
| Uniprot | Q92888 |
| GeneID | 9138; |

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.089ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ARHGEF1

Note: This product is for in vitro research use only