

Mouse Arhgef1/Rho guanine nucleotide exchange factor 1 ELISA Kit

Catalog No: #EK19755

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Arhgef1/Rho guanine nucleotide exchange factor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Rho guanine nucleotide exchange factor 1
Target Name	Rho guanine nucleotide exchange factor 1
Other Names	Lbc's second cousin,Lymphoid blast crisis-like 2,Lbcl2,Lsc
Accession No.	Q61210
Uniprot	Q61210
GeneID	16801;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:34pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Arhgef1

Note: This product is for in vitro research use only