

Rat Arhgef1/Rho guanine nucleotide exchange factor 1 ELISA Kit



Catalog No: #EK19756

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Arhgef1/Rho guanine nucleotide exchange factor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Rho guanine nucleotide exchange factor 1
Target Name	Rho guanine nucleotide exchange factor 1
Other Names	Lbc's second cousin,Lsc
Accession No.	Q9Z1I6
Uniprot	Q9Z1I6

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.089ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Arhgef1

Note: This product is for in vitro research use only