

Mouse Cfc1/Cryptic protein ELISA Kit

Catalog No: #EK19758

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cfc1/Cryptic protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cryptic protein
Target Name	Cryptic protein
Accession No.	P97766
Uniprot	P97766
GeneID	12627;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 20pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cfc1

Note: This product is for in vitro research use only