

Human ANKRD30A/Ankyrin repeat domain-containing protein 30A ELISA Kit



Catalog No: #EK19759

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ANKRD30A/Ankyrin repeat domain-containing protein 30A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Ankyrin repeat domain-containing protein 30A
Target Name	Ankyrin repeat domain-containing protein 30A
Other Names	Serologically defined breast cancer antigen NY-BR-1
Accession No.	Q9BXX3
Uniprot	Q9BXX3
GeneID	91074;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 41.4pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ANKRD30A

Note: This product is for in vitro research use only