

Mouse Adam33/Disintegrin and metalloproteinase domain-containing protein 33 ELISA Kit

Catalog No: #EK19760

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Adam33/Disintegrin and metalloproteinase domain-containing protein 33 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Disintegrin and metalloproteinase domain-containing protein 33
Target Name	Disintegrin and metalloproteinase domain-containing protein 33
Other Names	ADAM 33
Accession No.	Q923W9
Uniprot	Q923W9
GeneID	110751;

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:41.5pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Adam33

Note: This product is for in vitro research use only