

Rat Casp1/Caspase-1 ELISA Kit

Catalog No: #EK19763

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Casp1/Caspase-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Caspase-1
Target Name	Caspase-1
Other Names	CASP-1, Interleukin-1 beta convertase, IL-1BC, Interleukin-1 beta-converting enzyme, ICE, p45, Il1bc, Il1bce
Accession No.	P43527
Uniprot	P43527

Application Details

Assay Type: Sandwich Detect Range: 0.625-40ng/mL Sensitivity: 0.312ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Casp1

Note: This product is for in vitro research use only