

Human CSF1R/Macrophage colony-stimulating factor 1 receptor ELISA Kit



Catalog No: #EK19764

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CSF1R/Macrophage colony-stimulating factor 1 receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Macrophage colony-stimulating factor 1 receptor
Target Name	Macrophage colony-stimulating factor 1 receptor
Other Names	CSF-1 receptor,CSF-1-R,Proto-oncogene c-Fms,CD115,FMS
Accession No.	P07333
Uniprot	P07333
GeneID	1436;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CSF1R

Note: This product is for in vitro research use only