

Mouse Csf1r/Macrophage colony-stimulating factor 1 receptor ELISA Kit

Catalog No: #EK19765

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Csf1r/Macrophage colony-stimulating factor 1 receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Macrophage colony-stimulating factor 1 receptor
Target Name	Macrophage colony-stimulating factor 1 receptor
Other Names	CSF-1 receptor,CSF-1-R,Proto-oncogene c-Fms,CD115,Csfmr,Fms
Accession No.	P09581
Uniprot	P09581
GeneID	12978;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:38pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Csf1r

Note: This product is for in vitro research use only