

Rat Csf1r/Macrophage colony-stimulating factor 1 receptor ELISA Kit

Catalog No: #EK19766

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Csf1r/Macrophage colony-stimulating factor 1 receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Macrophage colony-stimulating factor 1 receptor
Target Name	Macrophage colony-stimulating factor 1 receptor
Other Names	CSF-1 receptor,CSF-1-R,Proto-oncogene c-Fms,CD115,Csfmr,Fms
Accession No.	Q00495
Uniprot	Q00495

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:41.6pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Csf1r

Note: This product is for in vitro research use only