

# Mouse Clec4e/C-type lectin domain family 4 member E ELISA Kit



Catalog No: #EK19767

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Mouse Clec4e/C-type lectin domain family 4 member E ELISA Kit                       |
| Brief Description  | ELISA Kits  |
| Applications       | ELISA   |
| Species Reactivity | Ms  |
| Specificity        | Natural and recombinant Mouse C-type lectin domain family 4 member E                |
| Target Name        | C-type lectin domain family 4 member E  |
| Other Names        | C-type lectin superfamily member 9,Macrophage-inducible C-type lectin,Clecsf9,Mincl |
| Accession No.      | Q9R0Q8  |
| Uniprot            | Q9R0Q8  |
| GeneID             | 56619;  |

## Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.08ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:Clec4e

Note: This product is for in vitro research use only