

Rat Clec4e/C-type lectin domain family 4 member E ELISA Kit

Catalog No: #EK19768

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Clec4e/C-type lectin domain family 4 member E ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat C-type lectin domain family 4 member E
Target Name	C-type lectin domain family 4 member E
Other Names	C-type lectin superfamily member 9,Macrophage-inducible C-type lectin,Clecsf9,Mincle
Accession No.	Q67EQ1
Uniprot	Q67EQ1
GeneID	450223;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.081ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Clec4e

Note: This product is for in vitro research use only